

**HACKETTSTOWN REGIONAL MEDICAL CENTER
CARDIO PULMONARY POLICY MANUAL
ROOM TEMPERATURE**

=====

Effective Date: July 2009

Cross Referenced:

Reviewed Date: 10/15

Revised Date: 06/12

Policy No: 4A.011

Origin: Cardio Pulmonary

Authority: Cardio/Pulmonary Manager

Page: 1 of 1

Purpose: The process in which to follow when recording the room temperature in accordance to manufacturer suggestions for cartridges.

Procedure: As part of the daily maintenance, record the room temperature to assure a safe temperature as recommended by Siemens Medical for the Wash/Waste cartridges. The recommended temperature range is between 2 degrees Celsius and 25 degrees Celsius. Record the temperature in the log. If the temperature is out of range, notify maintenance, so they can adjust the temperature.