

**HACKETTSTOWN REGIONAL MEDICAL CENTER  
CARDIO PULMONARY POLICY MANUAL  
DAILY CLEANING**

=====

**Effective Date: May 2007**

**Cross Referenced:**

**Reviewed Date: 02/16**

**Revised Date: 02/16**

**Policy No: 5.009**

**Origin: Cardio Pulmonary**

**Authority: Cardio/Pulmonary Manager**

**Page: 1 of 1**

---

**Purpose:** To assure that all patient related equipment is cleaned between each patient.

**Procedure:** All patient lead wires will be wiped down with Super Sani-cloth wipes between patients. If there is any adhesive residue on the wires, they must be cleaned with adhesive remover pads first, then wiped down with Super Sani-cloth wipes. Cables also must be wiped down as well. The Programmer is also wiped between patients.