

**HACKETTSTOWN REGIONAL MEDICAL CENTER
CARDIO PULMONARY POLICY MANUAL
DAILY CLEANING**

=====

Effective Date: March 2010

Cross Referenced:

Reviewed Date: 02/16

Revised Date: 12/12

Policy No: 8B.014

Origin: Cardio Pulmonary

Authority: Cardio/Pulmonary Manager

Page: 1 of 1

Purpose: To ensure that all patient equipment is cleaned between each patient.

Policy: All patient lead wires will be wiped down with Wet Task Wiping System between patients. If there is any adhesive residue on the wires, they must be cleaned with the adhesive remover pads first, then wiped down with Germicidal Wet Wipes (PDI-Sanicloth). Cables also must be wiped down as well. The treadmill is also wiped down between patients.