HACKETTSTOWN REGIONAL MEDICAL CENTER LABORATORY POLICY MANUAL CEREBRAL SPINAL FLUID COLLECTION REQUIREMENTS

Effective Date: March, 1998 Policy No: SC 100.09

Cross Referenced:

Reviewed Date: 5/2012

Origin: Special Collections
Authority: Laboratory Director

Revised Date: 1/12 Page: 1 of 3

PURPOSE: To describe conditions under which CSF is handled.

POLICY: As follows:

Departments: Chemistry and Hematology

Scheduling: Stat

Container: Sterile CSF lumbar puncture tubes, labeled 1, 2 and 3

Tests Included: Color

Appearance Xanthochromia Total WBC Count

Differential if the WBC is greater than or equal to 5/mm³

Total RBC Glucose Total Protein

Any other tests such as: Culture (Gram stain is automatically

performed with all CSF cultures), Bacterial Antigens, VDRL, Cytology, etc. must

be ordered separately.

Instructions:

- 1. Specimen must arrive in the laboratory within 30 minutes of collection. Deliver the specimen to specimen processing and bring the specimen to laboratory personnel's attention
- 2. All tubes are to be labeled with the patient's name, medical record number, doctor, admitting location, and date.
- 3. Tubes will be used in the following order unless specified by the physician on the tube, #1 for chemistry, #2 or #3 for cell count and culture.
- 4. Stat gram stain performed and results called to doctor.
- 5. Culture and/or Bacterial Antigens sent stat to St. Clare's. Record results of gram stain on requisition.