Section: Division of Nursing

\*\*\*\*\*\*\*\*\*\*\*\*
\* PROCEDURE \*

Index: 6080.007a Page: 1 of 1

Issue Date: February 26, 1992 Review Date: July 23, 2010

HACKETTSTOWN REGIONAL MEDICAL CENTER

Originator: K. Shubick, RN

3-NORTH (Scope)

TITLE: APPLICATION OF CERVICAL TRACTION PROCEDURE

**PURPOSE:** To outline the process for applying cervical traction.

**EQUIPMENT:** 1. Cervical traction bracket

Head halter
 Spreader bar
 Weight carrier
 Traction cord
 Weights

CONTENT: PROCEDURE KEY POINTS

 Attach cervical traction bracket to foot at head of bed.

2. Turn bed around so that foot is at the head.

- 3. Tie traction cord to spreader bar and feed through pulleys and tie to weight carrier.
- 4. Apply head halter to patient.
- 5. Attach head halter to spreader bar.
- 6. Apply weights.
- Document type of traction applied, amount of weight used, and patient's response to traction in nursing notes.