

Section: Division of Nursing

* **PROCEDURE** *

Index: 7420.078a
Page: 1 of 1
Issue Date: May 23, 2005

Approval: _____

HACKETTSTOWN REGIONAL MEDICAL CENTER

Originator: D. Greco, RN
Reviewed by: W. Zeiher, RN

OR
(Scope)

TITLE: ESSURE MICRO INSERTION

PURPOSE: To provide instructions for the Essure Micro Insertion Procedure.

EQUIPMENT LIST:

1. Hot line machine & hot line tubing
2. I/V extension set X2
3. Essure Micro Insertion Kit
(Found in GYN section – Central Supply)
4. I/V pressure bag
5. Storz tower

CONTENT: PROCEDURE STEPS:

1. Plug hot line pump in and prime hot line tubing in order to insure warming prior to physician and patient entering room.
2. Attach pressure bag to I/V 1000ML normal saline and pump bag to 80mmHg (*must be maintained @ 80mmHg – no higher unless specifically ordered by physician).
3. Place patient on bed and place in lithotomy using candy cane stirrups to secure legs.
4. Place Storz Tower #4 on patient's right side near legs.
5. Extension tubing will be passed by surgeon to circulator which attached to hot line tubing containing warmed fluid (Sterile end attached to hysteroscope).
6. Open smaller package from Essure Kit containing the dilator and place on sterile field using aseptic technique.
7. Using aseptic technique open one long package containing the Essure delivery system and hand to scrub nurse.
8. After physician completes insertion open second long packet and hand to scrub nurse in same aseptic manner.
9. Make sure that picture is taken of each insertion and given to surgeon.
10. Complete patient identification card and include a patient information booklet.

Reference: Conceptus Incorporated Package Insert, 2004.