

Section: Division of Nursing
Approval: _____

PROCEDURE

Index: 7420.071a
Page: 1 of 1
Issue Date: May 30, 2002
Review Date: March 12, 2010

HACKETTSTOWN REGIONAL MEDICAL CENTER

Originator: M. Mason, RN
Revised by: M. Mason, RN
Reviewed by: R. Trimmer, RN

OR
(Scope)

TITLE: SENS Rx GAMMA FINDER

PURPOSE: To outline the procedure for set-up and use of Cath-Finder Catheter Tracking System.

SUPPORTIVE DATA: The Deltic Cath-Finder Catheter Tracking System is designed for placing and locating Deltic implantable and/or pcutaneous intravenous Catheters during the Catheter placement procedure.

EQUIPMENT LIST:

1. One locator wand
2. Sterile Sleeve

CONTENT

PROCEDURE STEPS

KEY POINTS

1. Open sterile cover on bluid.
2. Put locator wand in sterile cover.
3. Turn the locator wand on by pressing the locator wand's button. Catheter/sensor assembly to verify operation.

Cleaning:

Clean wand with damp soft cloth and mild detergent.

Warning:

Do not autoclave, irradiate, or gas-sterilize the electronic controller, locator wand, or extension cables.

References: Seno Rx, Inc., 11 Columbia, Suite A, Akso Viejo, CA 92656